



Guarantee high purity and yield from virus sample including 'Novel H1N1 Flu'

Viral DNA/RNA extraction ability equivalent to supplier Q


Feline herpesvirus (clinical sample)




Canine distemper virus (vaccine)



Newcastle disease virus (clinical sample)



Application for viral RNA isolation of 'Novel H1N1 Flu'



Result of 'Novel H1N1 Flu' detection with viral RNA from clinical sample using **INTRON Viral Gene-spin™ kit** and 'Novel H1N1 Flu' detection using company S's FluA subtyping kit.


Clean Automatic Packaging System (CAPS)




Extract viral RNA from almost virus samples

Serum, Plasma, Whole blood, Tissue, Swab, Urine, Cultured cell, Cell-free fluid, Virus infection tissue or cell and etc..


The test results of viral RNA extraction from a various of virus using Viral Gene-spin™ Kit



Fast extraction compared to a competitor by simple process



Simple and Easy Protocol



Reference

- Characterization of a Recombinant Newcastle Disease Virus Vaccine Strain *Ann. Clin. Lab. Sci.*, (2008) 38: 41 - 46.
- Comparison of the Seeplex Reverse Transcription PCR Assay with the R-mix Viral Culture and Immunofluorescence Techniques for Detection of Eight Respiratory Viruses *Avian Dis.* (2007) 51 (2):540-5
- Development of one-step reverse transcriptase-polymerase chain reaction to detect duck hepatitis virus type 1. *J. Virol Methods* (2009) 156(1-2) : 111-116
- Rapid detection and identification of 12 respiratory viruses using a dual priming oligonucleotide system-based multiplex PCR assay *J. Vet. Sci.* (2008) 9(2), 207-209
- Prevalence of feline herpesvirus 1, feline calicivirus and Chlamydomphila felis in clinically normal cats at a Korean animal shelter *J. Vet. Sci.* (2009) 10(2), 131-139
- Protection of chicken against very virulent IBDV provided by in ovo priming with DNA vaccine and boosting with killed vaccine and the adjuvant effects of plasmid-encoded chicken interleukin-2 and interferon- γ *J. Microbiology* (2008) 709-719
- Molecular characterization of ORFs 2 to 7 of Korean porcine reproductive and respiratory syndrome virus (CA) and its protein expression by recombinant baculoviruses *Laryngoscope.* (2007) 117(8):1445-9
- Detection of Parainfluenza Virus 3 in Turbinate Epithelial Cells of Postviral Olfactory Dysfunction Patients *J. Virol Methods* (2008) 149 : 76-84
- Direct detection of lamivudine-resistant hepatitis B virus mutants by a multiplex PCR using dual-priming oligonucleotide primers *Virus Genes* (2007) 35:339-345
- Prevalence and molecular characterization of meq in feather follicular epithelial cells of Korean broiler chickens Detection of rhinovirus in turbinate epithelial cells of chronic sinusitis *Am J Rhinol* (2006) 20 : 634-636
- *OIE Terrestrial Manual* (2008) Chapter 2. 3. 8. Duck virus hepatitis